

RABEPK Polyclonal Antibody

Catalog No. E-AB-19004

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Fusion protein of human RABEPK
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4


Applications Recommended Dilution

IHC	1:50-1:300
------------	------------

Data



Immunohistochemistry of paraffin-embedded Human tonsil tissue using RABEPK Polyclonal Antibody at dilution of 1:85(×200)



Immunohistochemistry of paraffin-embedded Human brain tissue using RABEPK Polyclonal Antibody at dilution of 1:85(×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Rab9 effector required for endosome to trans-Golgi network (TGN) transport.

For Research Use Only