

GPLOW Polyclonal Antibody

Catalog No. E-AB-19058

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Fusion protein of human GPLOW
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution


WB	1:1000-1:5000
IHC	1:100-1:500

Data



Western blot analysis of 293T and Jurkat cell lysates using GPLOW Polyclonal Antibody at dilution of 1:2000

Observed Mw: Refer to figures
Calculated Mw: 53 kDa



Immunohistochemistry of paraffin-embedded Human tonsil tissue using GPLOW Polyclonal Antibody at dilution of 1:140 (x200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

This gene encodes a putative RNA-binding protein containing G-patch and KOW (Kyprides, Ouzounis, Woese) domains. The encoded protein interacts directly with protein kinase A and protein kinase X and is also found associated with the spliceosome.

For Research Use Only