

SARS Polyclonal Antibody

Catalog No. E-AB-19068

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Fusion protein of human SARS
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol,pH7.4

Applications Recommended Dilution


WB	1:1000-1:5000
IHC	1:50-1:300

Data



Western blot analysis of Human fetal brain tissue lysate using SARS Polyclonal Antibody at dilution of 1:850

Observed Mw:Refer to figures
Calculated Mw:59 kDa



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using SARS Polyclonal Antibody at dilution of 1:65(×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

This gene belongs to the class II amino-acyl tRNA family. The encoded enzyme catalyzes the transfer of L-serine to tRNA (Ser) and is related to bacterial and yeast counterparts. Multiple alternatively spliced transcript variants have been described but the biological validity of all variants is unknown.

For Research Use Only