

DPY19L3 Polyclonal Antibody

Catalog No. E-AB-19096

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Fusion protein of human DPY19L3
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol,pH7.4

Applications Recommended Dilution

IHC	1:100-1:300
------------	-------------

Data



Immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using DPY19L3 Polyclonal Antibody at dilution of 1:120(×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Probable C-mannosyltransferase that mediates C-mannosylation of tryptophan residues on target proteins.

For Research Use Only