

PHLDA3 Polyclonal Antibody

Catalog Number:E-AB-19111

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Fusion protein of human PHLDA3
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol,pH7.4


Applications Recommended Dilution

IHC	1:100-1:300
ELISA	1:5000-1:10000

Data



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using PHLDA3 Polyclonal Antibody at dilution of 1:110(×200)



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using PHLDA3 Polyclonal Antibody at dilution of 1:110(×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

PHLDA3 (Pleckstrin Homology Like Domain Family A Member 3) is a Protein Coding gene. Among its related pathways are Apoptosis and Autophagy and PI3K / Akt Signaling. GO annotations related to this gene include phosphatidylinositol-4,5-bisphosphate binding and phosphatidylinositol-3-phosphate binding. An important paralog of this gene is PHLDA1.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017