

IRF2 Polyclonal Antibody

Catalog No. E-AB-19511

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Synthetic peptide of human IRF2
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution

WB 1:1000-1:5000

Data



Western blot analysis of Raji and Jurkat cell lysates using IRF2 Polyclonal Antibody at dilution of 1:1350

Observed Mw: Refer to figures

Calculated Mw: 39 kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

IRF2 encodes interferon regulatory factor 2, a member of the interferon regulatory transcription factor (IRF) family. IRF2 competitively inhibits the IRF1-mediated transcriptional activation of interferons alpha and beta, and presumably other genes that employ IRF1 for transcription activation. However, IRF2 also functions as a transcriptional activator of histone H4.

For Research Use Only