

NIPSNAP2 Polyclonal Antibody

Catalog No. E-AB-19832

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Synthetic peptide of human NIPSNAP2
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution

WB	1:500-1:2000
-----------	--------------

Data



Western blot analysis of HeLa HT-29 and HEPG2 cell
Human liver tissue Human fetal liver tissue Mouse
liver tissue and Human fetal brain tissue lysates using
NIPSNAP2 Polyclonal Antibody at dilution of 1:400

Observed Mw: Refer to figures

Calculated Mw: 34 kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

This gene encodes a member of the NipSnap family of proteins that may be involved in vesicular transport. The encoded protein is localized to mitochondria and plays a role in oxidative phosphorylation. A pseudogene of this gene is located on the long arm of chromosome 2. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

For Research Use Only