

## USP50 Polyclonal Antibody

**Catalog No.** E-AB-52022

*Note:* Centrifuge before opening to ensure complete recovery of vial contents.


### Description

<b>Reactivity</b>	Human, Mouse
<b>Immunogen</b>	Synthetic peptide of human USP50
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification</b>	Antigen affinity purification
<b>Conjugation</b>	Unconjugated
<b>Buffer</b>	PBS with 0.05% NaN <sub>3</sub> and 40% Glycerol,pH7.4

### Applications Recommended Dilution

**WB** 1:500-1:2000

### Data



Western blot analysis of Mouse adrenal gland tissue lysate using USP50 Polyclonal Antibody at dilution of 1:400

**Observed Mw:Refer to figures**  
**Calculated Mw:39 kDa**

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### Background

USP50 (Ubiquitin Specific Peptidase 50) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include thiol-dependent ubiquitinyl hydrolase activity. An important paralog of this gene is USP8.

### For Research Use Only