

Phospho-MST1R (Ser1394) Polyclonal Antibody

Catalog No. E-AB-21033

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Synthesized peptide derived from human Ron around the phosphorylation site of S1394.
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Buffer	PBS with 0.02% sodium azide, 0.5% protective protein and 50% glycerol, pH7.4

Applications Recommended Dilution

WB	1:500-1:2000
ELISA	1:40000

Data



Western Blot analysis of 3T3 cells using Phospho-MST1R (Ser1394) Polyclonal Antibody at dilution of 1:500.

Observed Mw:152kDa
Calculated Mw:152kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Receptor for macrophage stimulating protein (MSP). Has a tyrosine-protein kinase activity.

For Research Use Only