

Phospho-PDGFR α (Tyr849) Polyclonal Antibody

Catalog No. E-AB-21070

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Synthesized peptide derived from human PDGFR- α around the phosphorylation site of Y849.
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Buffer	PBS with 0.02% sodium azide, 0.5% protective protein and 50% glycerol, pH7.4

Applications Recommended Dilution

WB	1:500-1:2000
IHC	1:100-1:300

Data



Western Blot analysis of 293 cells with Phospho-PDGFR α (Tyr849) Polyclonal Antibody
Observed Mw:140kDa
Calculated Mw:123kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Receptor that binds both PDGFA and PDGFB and has a tyrosine-protein kinase activity.

For Research Use Only