

# Strep-Tag Monoclonal Antibody

Catalog Number:E-AB-20021



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

## Description

<b>Immunogen</b>	Synthetic Peptide
<b>Host</b>	Mouse
<b>Isotype</b>	IgG
<b>Clone</b>	Clone:2L3
<b>Purification</b>	Protein A purification
<b>Conjugation</b>	Unconjugated
<b>Formulation</b>	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

## Applications Recommended Dilution

<b>WB</b>	1:5000-1:10000
-----------	----------------

## Data



www.elabscience.com

Western Blot analysis of 1ug Strep II fusion protein using Strep-Tag Monoclonal Antibody at dilution of 1) 1:5000 2) 1:10000.

## Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017