

S-Tag Monoclonal Antibody

Catalog No. E-AB-20022

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Immunogen	Synthetic Peptide
Host	Mouse
Isotype	IgG
Clone	Clone:8C3
Purification	Protein A purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB	1:5000-1:10000
-----------	----------------

Data



www.elabscience.com

Western Blot analysis of 2ug S-Tag fusion protein using S-Tag Monoclonal Antibody at dilution of 1) 1:5000 2) 1:10000.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only