

GST-Tag Monoclonal Antibody

Catalog No. E-AB-20075

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Immunogen	Recombinant Protein
Host	Mouse
Isotype	IgG
Clone	Clone:9A3
Purification	Protein A purification
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB	1:2000-5000
-----------	-------------

Data



Western Blot analysis of 0.2ug GST protein, using GST Monoclonal Antibody at dilution of 1:5000.

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only