

RELB Polyclonal Antibody

Catalog Number:E-AB-40045

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Rat
Immunogen	Recombinant human Transcription factor RelB protein
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Formulation	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.


Applications Recommended Dilution

IHC	1:50-1:100
------------	------------

Data



Immunohistochemistry of paraffin-embedded Rat spleen using RELB Polyclonal Antibody at dilution of 1:50



Immunohistochemistry analysis of paraffin-embedded Human Kidney using RELB Polyclonal Antibody(Elabscience® Product Detected by Lifespan).

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles.
----------------	---

Background

RELB,also named as Transcription factor RelB,is a 579 amino acid protein, which contains one RHD domain.RELB localizes in the nucleus and cytoplasm. RELB is a pleiotropic transcription factor which is present in almost all cell types and is involved in many biological processes such as inflammation, immunity,differentiation, cell growth, tumorigenesis and apoptosis.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017