

# ANXA5 Polyclonal Antibody

Catalog Number:E-AB-40114

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


## Description

<b>Reactivity</b>	Human,Mouse,Rat
<b>Immunogen</b>	Recombinant human Annexin A5 protein
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification</b>	Affinity purification
<b>Conjugation</b>	Unconjugated
<b>Formulation</b>	PBS with 0.05% Proclin300 and 50% glycerol, pH7.4.

## Applications Recommended Dilution


<b>WB</b>	1:500-1:1000
<b>IHC</b>	1:50-1:100

## Data




Western Blot analysis of HeLa and Mouse kidney using ANXA5 Polyclonal Antibody at dilution of 1:600

**Observed Mw:36KDa**  
**Calculated Mw:36KDa**



Immunohistochemistry of paraffin-embedded Rat kidney using ANXA5 Polyclonal Antibody at dilution of 1:50



Immunohistochemistry analysis of paraffin-embedded Human Prostate using ANXA5 Polyclonal Antibody(Elabscience® Product Detected by Lifespan).

## Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

## Background

This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017

# ANXA5 Polyclonal Antibody

Catalog Number:E-AB-40114



involved in the blood coagulation cascade.

---

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017