

TRAPPC3 Polyclonal Antibody

Catalog No. E-AB-30647

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Synthesized peptide derived from the C-terminal region of human BET3
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Buffer	PBS with 0.02% sodium azide,0.5% protective protein and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB	1:500-1:2000
IHC	1:100-1:300
ELISA	1:5000

Data



Western Blot analysis of COLO205 cells using
TRAPPC3 Polyclonal Antibody at dilution of 1:500.

Observed Mw:22kDa
Calculated Mw:20kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

May play a role in vesicular transport from endoplasmic reticulum to Golgi.

For Research Use Only