

SHC3 Polyclonal Antibody

Catalog No. E-AB-52214

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Fusion protein of human SHC3
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4


Applications Recommended Dilution

WB	1:500-1:2000
IHC	1:50-1:100

Data



Western blot analysis of PC3 cell lysate using SHC3 Polyclonal Antibody at dilution of 1:650
Observed Mw: Refer to figures
Calculated Mw: 64 kDa



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using SHC3 Polyclonal Antibody at dilution of 1:40 (×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Signaling adapter that couples activated growth factor receptors to signaling pathway in neurons. Involved in the signal transduction pathways of neurotrophin-activated Trk receptors in cortical neurons.

For Research Use Only