

GATD3B Polyclonal Antibody

Catalog No. E-AB-52500

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Fusion protein of human GATD3B
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol,pH7.4

Applications Recommended Dilution

WB 1:500-1:2000

Data



Western blot analysis of 293T cell lysates using GATD3B Polyclonal Antibody at dilution of 1:250

Observed Mw:Refer to figures
Calculated Mw:28 kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

GATD3B (Glutamine Amidotransferase Like Class 1 Domain Containing 3B) is a Protein Coding gene. Among its related pathways are Gastric Cancer Network 1. An important paralog of this gene is GATD3A.

For Research Use Only