

EIF2S3 Polyclonal Antibody

Catalog No. E-AB-52647

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Rat
Immunogen	Fusion protein of human EIF2S3
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4


Applications Recommended Dilution

WB 1:500-1:2000,
IHC 1:30-1:150,
ELISA
1:5000-1:10000

Data



Western blot analysis of HL-60 cell lysate using
EIF2S3 Polyclonal Antibody at dilution of 1:450
Observed Mw:Refer to figures
Calculated Mw:51 kDa



Immunohistochemistry of paraffin-embedded
Human cervical cancer tissue using EIF2S3
Polyclonal Antibody at dilution of 1:60(×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

The protein encoded by this gene is the largest subunit of a heterotrimeric GTP-binding protein involved in the recruitment of methionyl-tRNA(i) to the 40 S ribosomal subunit.

For Research Use Only