

ACSF2 Polyclonal Antibody

Catalog No. E-AB-52654

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Fusion protein of human ACSF2
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4


Applications Recommended Dilution

WB	1:500-1:2000
IHC	1:30-1:150

Data



Western blot analysis of 293T cell Rat liver tissue A172 and Jurkat cell lysates using ACSF2 Polyclonal Antibody at dilution of 1:400
Observed Mw: Refer to figures
Calculated Mw: 68 kDa



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ACSF2 Polyclonal Antibody at dilution of 1:30 (×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Acyl-CoA synthases catalyze the initial reaction in fatty acid metabolism, by forming a thioester with CoA. Has some preference toward medium-chain substrates. Plays a role in adipocyte differentiation.

For Research Use Only