

## CCDC4 Polyclonal Antibody

Catalog No. E-AB-30792

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


### Description

<b>Reactivity</b>	Human
<b>Immunogen</b>	Synthesized peptide derived from the Internal region of human CCDC4
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Purification</b>	Affinity purification
<b>Buffer</b>	PBS with 0.02% sodium azide, 0.5% protective protein and 50% glycerol pH 7.4.

### Applications Recommended Dilution

<b>WB</b>	1:500-1:2000
<b>ELISA</b>	1:40000

### Data



Western Blot analysis of HuvEc cells with CCDC4  
Polyclonal Antibody.  
**Observed Mw:57kDa**  
**Calculated Mw:58kDa**

### Preparation & Storage

**Storage** Store at -20°C. Avoid freeze / thaw cycles.

### Background

BEND4 (BEN Domain Containing 4) is a Protein Coding gene.

### For Research Use Only