

PDAP1 Polyclonal Antibody

Catalog No. E-AB-52832

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Fusion protein of human PDAP1
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution

WB	1:500-1:2000
IHC	1:50-1:300

Data



Western blot analysis of 293T cell and Rat brain tissue lysates using PDAP1 Polyclonal Antibody at dilution of 1:400

Observed Mw: Refer to figures
Calculated Mw: 21 kDa



Immunohistochemistry of paraffin-embedded Human tonsil tissue using PDAP1 Polyclonal Antibody at dilution of 1:50 (×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

The protein encoded by this gene is a phosphoprotein that may upregulate the PDGFA-stimulated growth of fibroblasts and also downregulate the mitogenicity of PDGFB. The encoded protein in rodents has been shown to bind PDGFA with a low affinity.

For Research Use Only