

AKAP17A Polyclonal Antibody

Catalog No. E-AB-52956

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Fusion protein of human AKAP17A
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol,pH7.4

Applications Recommended Dilution

IHC	1:40-1:200
------------	------------

Data



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using AKAP17A Polyclonal Antibody at dilution of 1:80(×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

This locus encodes a protein kinase A anchoring protein. The encoded protein is part of the spliceosome complex and is involved in the regulation of alternate splicing in some mRNA precursors. Alternatively spliced transcript variants have been identified for this gene. AKAP17A (A-Kinase Anchoring Protein 17A) is a Protein Coding gene. Diseases associated with AKAP17A include Null-Cell Leukemia and Chronic Tic Disorder. GO annotations related to this gene include poly(A) RNA binding and protein kinase A binding.

For Research Use Only