

FDXR Polyclonal Antibody

Catalog No. E-AB-52975

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Fusion protein of human FDXR
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution


WB	1:1000-1:5000
IHC	1:100-1:500

Data




Western blot analysis of HepG2 cell lysate using FDXR Polyclonal Antibody at dilution of 1:1000

Observed Mw: Refer to figures
Calculated Mw: 54 kDa



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using FDXR Polyclonal Antibody at dilution of 1:110 (×200)



Immunohistochemistry of paraffin-embedded Human tonsil tissue using FDXR Polyclonal Antibody at dilution of 1:110 (×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

For Research Use Only

Background

This gene encodes a mitochondrial flavoprotein that initiates electron transport for cytochromes P450 receiving electrons from NADPH. Multiple alternatively spliced transcript variants have been found for this gene.

For Research Use Only