

MYEF2 Polyclonal Antibody

Catalog No. E-AB-53086

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Fusion protein of human MYEF2
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution


WB	1:500-1:2000
IHC	1:50-1:200

Data




Western blot analysis of A172 and TM4 cell lysates using MYEF2 Polyclonal Antibody at dilution of 1:800

Observed Mw: Refer to figures
Calculated Mw: 64 kDa



Immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using MYEF2 Polyclonal Antibody at dilution of 1:70 (×200)



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using MYEF2 Polyclonal Antibody at dilution of 1:70 (×200)

Preparation & Storage

For Research Use Only

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Transcriptional repressor of the myelin basic protein gene (MBP). Binds to the proximal MB1 element 5'-TTGTCC-3' of the MBP promoter. Its binding to MB1 and function are inhibited by PURA (By similarity).

For Research Use Only