

SUV39H2 Polyclonal Antibody

Catalog No. E-AB-53163

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse
Immunogen	Fusion protein of human SUV39H2
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution

IHC	1:50-1:200
------------	------------

Data



Immunohistochemistry of paraffin-embedded Human prostate cancer tissue using SUV39H2 Polyclonal Antibody at dilution of 1:70 (×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

SUV39H2 (Suppressor Of Variegation 3-9 Homolog 2) is a Protein Coding gene. Among its related pathways are Circadian rhythm related genes and Lysine degradation. GO annotations related to this gene include chromatin binding and transcription regulatory region sequence-specific DNA binding. An important paralog of this gene is SUV39H1.

For Research Use Only