

CD5 Polyclonal Antibody

Catalog No. E-AB-30831

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Synthesized peptide derived from human CD5 around the non-phosphorylation site of Y453.
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Buffer	PBS with 0.02% sodium azide,0.5% protective protein and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB	1:500-1:2000
ELISA	1:20000

Data



Western Blot analysis of 293 cells with CD5 Polyclonal Antibody.

Observed Mw:60kDa
Calculated Mw:55kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2.

For Research Use Only