

BRSK2 Polyclonal Antibody

Catalog No. E-AB-53606

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Rat
Immunogen	Synthetic peptide of human BRSK2
Host	Rabbit
Isotype	IgG
Purification	Antigen affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4

Applications Recommended Dilution

IHC	1:40-1:200
------------	------------

Data



Immunohistochemistry of paraffin-embedded Human brain tissue using BRSK2 Polyclonal Antibody at dilution of 1:70(×200)

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

BRSK2 (BR Serine/Threonine Kinase 2) is a Protein Coding gene. Diseases associated with BRSK2 include Limbic Encephalitis and Sporotrichosis. Among its related pathways are DNA Damage and LKB1 signaling events. Gene Ontology (GO) annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is BRSK1

For Research Use Only