

CERKL Polyclonal Antibody

Catalog No. E-AB-30890

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Synthesized peptide derived from the Internal region of human CERKL
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Buffer	PBS with 0.02% sodium azide, 0.5% protective protein and 50% glycerol, pH7.4

Applications Recommended Dilution

WB	1:500-1:2000
IHC	1:100-1:300
ELISA	1:20000

Data



Western Blot analysis of A549 cells with CERKL
Polyclonal Antibody.
Observed Mw:63kDa
Calculated Mw:63kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Has no detectable ceramide-kinase activity. Overexpression of CERKL protects cells from apoptosis in oxidative stress conditions.

For Research Use Only