

CCZ1/CCZ1B Polyclonal Antibody

Catalog No. E-AB-30903

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human, Mouse, Monkey
Immunogen	Synthesized peptide derived from the C-terminal region of human CGI-43
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Buffer	PBS with 0.02% sodium azide, 0.5% protective protein and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB	1:500-1:2000
ELISA	1:10000

Data



Western Blot analysis of Jurkat cells using
CCZ1/CCZ1B Polyclonal Antibody at dilution of
1:500.

Observed Mw:56kDa
Calculated Mw:56kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

CCZ1 (CCZ1 Homolog, Vacuolar Protein Trafficking And Biogenesis Associated) is a Protein Coding gene. Among its related pathways are Vesicle-mediated transport and RAB GEFs exchange GTP for GDP on RABs. An important paralog of this gene is CCZ1B.

For Research Use Only