

CLCC1 Polyclonal Antibody

Catalog No. E-AB-30951

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Synthesized peptide derived from the Internal region of human CLCC1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.02% sodium azide, 0.5% protective protein and 50% glycerol, pH7.4

Applications Recommended Dilution

WB	1:500-1:2000
IF	1:100-1:300
ELISA	1:40000

Data



Western Blot analysis of K562, Jurkat cells using CLCC1 Polyclonal Antibody at dilution of 1:1000.

Observed Mw:62kDa
Calculated Mw:62kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

CLCC1 (Chloride Channel CLIC Like 1) is a Protein Coding gene. GO annotations related to this gene include chloride channel activity.

For Research Use Only