

CLIP4 Polyclonal Antibody

Catalog No. E-AB-30957

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Synthesized peptide derived from the Internal region of human CLIP4
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Buffer	PBS with 0.02% sodium azide, 0.5% protective protein and 50% glycerol pH 7.4.

Applications Recommended Dilution

WB	1:500-1:2000
ELISA	1:40000

Data



Western Blot analysis of HepG2 cells with CLIP4
Polyclonal Antibody.
Observed Mw:76kDa
Calculated Mw:76kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

CLIP4 (CAP-Gly Domain Containing Linker Protein Family Member 4) is a Protein Coding gene. An important paralog of this gene is CLIP3.

For Research Use Only