

CD272 Polyclonal Antibody

Catalog No. E-AB-11056

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description


Reactivity	Human
Immunogen	Recombinant protein of human BTLA
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% sodium azide and 50% glycerol, PH7.4

Applications


Recommended Dilution

IHC 1:25-1:100


Data



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using CD272 Polyclonal Antibody at dilution 1:30



Immunohistochemistry of paraffin-embedded Human tonsil tissue using CD272 Polyclonal Antibody at dilution 1:30



Immunohistochemistry of paraffin-embedded Small Intestine tissue using CD272 Polyclonal Antibody at dilution of 1:50 (Elabscience® Product Detected by Lifespan).

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

For Research Use Only

This gene encodes a member of the immunoglobulin superfamily. The encoded protein contains a single immunoglobulin (Ig) domain and is a receptor that relays inhibitory signals to suppress the immune response. Alternative splicing results in multiple transcript variants. Polymorphisms in this gene have been associated with an increased risk of rheumatoid arthritis.