

EMC8 Polyclonal Antibody

Catalog No. E-AB-14013

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant protein of human EMC8
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% sodium azide and 50% glycerol, PH7.4


Applications Recommended Dilution

WB	1:500-1:2000
IHC	1:50-1:200


Data



Western Blot analysis of HeLa and HT-29 cell using EMC8 Polyclonal Antibody at dilution of 1:650
Calculated Mw:24kDa



Immunohistochemistry of paraffin-embedded Human cervical cancer using EMC8 Polyclonal Antibody at dilution of 1:30



Immunohistochemistry of paraffin-embedded Human brain using EMC8 Polyclonal Antibody at dilution of 1:30

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

For Research Use Only

The protein encoded by this gene has ubiquitin ligase activity. It mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including TP53, HDAC1 and CDKN1B, thus regulating their levels and cell cycle progression. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.