

HMGB4 Polyclonal Antibody

Catalog No. E-AB-14128

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Recombinant protein of human HMGB4
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% sodium azide and 50% glycerol, PH7.4

Applications Recommended Dilution

WB 1:200-1:1000

Data



Western Blot analysis of 231 cell and Human fetal liver tissue using HMGB4 Polyclonal Antibody at dilution of 1:500

Calculated Mw:22kDa

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

High mobility group protein B4 is a transcription factor that in humans is encoded by the HMGB4 gene. HMGB4 is strongly and preferentially expressed in the adult mouse testis and weakly in the brain, but not in many other tissues. HMGB4 associates with chromatin, and in transfection assays, in contrast to HMGB1, it acts as a potent transcriptional repressor. During spermatogenesis, HMGB4 is present in the euchromatin of late pachytene spermatocytes and haploid round spermatids, whereas stronger expression is observed during the elongation phase, where it localizes to the basal pole of the nucleus in a manner mutually exclusive with H1FNT (H1T2) localized at the apical pole. HMGB4 basal localization is lost in H1FNT-mutant spermatids, showing that H1FNT provides a positional cue for organizing chromatin domains within the nucleus. These results show that HMGB4 and H1FNT specify distinct chromatin domains at the apical and basal poles of the elongating spermatid nucleus.

For Research Use Only