

C1QTNF6 Polyclonal Antibody

Catalog No. E-AB-14425

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Recombinant protein of human C1QTNF6
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% sodium azide and 50% glycerol, PH7.4

Applications Recommended Dilution

IHC 1:25-1:100

Data



Immunohistochemistry of paraffin-embedded Human brain tissue using C1QTNF6 Polyclonal Antibody at dilution 1:30

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

This protein is not uniformly distributed but rather focally localized in hepatocellular carcinoma cells and the tumor endothelial cells. It may play a role in metabolism and immunity. CTRPs acts as new molecular mediators connecting inflammatory and metabolic diseases.

For Research Use Only