

HEXIM1 Polyclonal Antibody

Catalog No. E-AB-14432

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human,Mouse,Rat
Immunogen	Recombinant protein of human HEXIM1
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% sodium azide and 50% glycerol, PH7.4


Applications Recommended Dilution

IHC 1:100-1:300

Data



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using HEXIM1 Polyclonal Antibody at dilution 1:65



Immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using HEXIM1 Polyclonal Antibody at dilution 1:65

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Expression of this gene is induced by hexamethylene-bis-acetamide in vascular smooth muscle cells. This gene has no introns.

For Research Use Only