

AQP11 Polyclonal Antibody

Catalog No. E-AB-17682

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Reactivity	Human
Immunogen	Synthetic peptide of human AQP11
Host	Rabbit
Isotype	IgG
Purification	Affinity purification
Conjugation	Unconjugated
Buffer	PBS with 0.05% sodium azide and 50% glycerol, PH7.4

Applications Recommended Dilution

IHC	1:20-100
ELISA	1:2000-5000

Data



Immunohistochemistry of paraffin-embedded Human liver cancer tissue using AQP11 Polyclonal Antibody at dilution 1:30

Preparation & Storage

Storage Store at -20°C. Avoid freeze / thaw cycles.

Background

Aquaporins facilitate the transport of water and small neutral solutes across cell membranes.

For Research Use Only