

Recombinant Human ERN1/IRE1 Protein (aa 465-977, His & GST Tag)

Catalog No. PKSH030346

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	hIRE1p;IRE1;IRE1a;IRE1P
Species	Human
Expression Host	Baculovirus-Insect Cells
Sequence	Pro 465-Leu 977
Accession	O75460-1
Calculated Molecular Weight	86.0 kDa
Observed molecular weight	90 kDa
Tag	N-His-GST
Bioactivity	The specific activity was determined to be 5 nmol/min/mg using MBP as substrate

Properties

Purity	> 80 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Shipping	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < -20°C.
Formulation	Supplied as sterile solution of 20mM Tris, 500mM NaCl, 10% glycerol, pH 7.4
Reconstitution	Not Applicable

Data



> 80 % as determined by reducing SDS-PAGE.

Background

The encoded protein contains two functional catalytic domains; a serine/threonine-protein kinase domain and an endoribonuclease domain. This protein functions as a sensor of unfolded proteins in the endoplasmic reticulum (ER) and triggers an intracellular signaling pathway termed the unfolded protein response (UPR). The UPR is an ER stress response that is conserved from yeast to mammals and activates genes involved in degrading misfolded proteins; regulating protein synthesis and activating molecular chaperones. This protein specifically mediates the splicing and activation of the stress response transcription factor X-box binding protein 1.

For Research Use Only