

Recombinant Human FABP3 Protein

Catalog Number:PKSH030830



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Fatty Acid-Binding Protein Heart;Fatty Acid-Binding Protein 3;Heart-Type Fatty Acid-Binding Protein;H-FABP;Mammary-Derived Growth Inhibitor;MDGIMuscle Fatty Acid-Binding Protein;M-FABP;FABP3;FABP11;MDGI;H-FABP;O-FABP
Species	Human
Expression Host	E.coli
Sequence	Met 1-Ala 133
Accession	P05413
Calculated Molecular Weight	14.9 kDa
Tag	None

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	Please contact us for more information.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile 50mM Tris, pH 8.0 Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Fatty acid binding protein 3 (FABP3; also termed heart-type fatty acid binding protein) is a member of the intracellular lipid-binding protein family that may be essential in fatty acid transport; cell growth; cellular signaling and gene transcription. Previously FABP3 was involved in apoptosis-associated congenital cardiac malformations. FABP3 knockdown exhibited significant toxic effects on cardiac development and mitochondrial function; which may be responsible for the knockdown of FABP3-induced apoptosis. FABP3 as a candidate gene underlying the etiology of congenital heart defects. Overexpression of FABP3 inhibited cell growth and proliferation via negative regulation of the cell cycle and down-regulation of cell growth factors; but enhances cell survival in hypoxic or ischemic conditions. FABPs are known to be carrier proteins for transporting fatty acids and other lipophilic substances from the cytoplasm to the

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Human FABP3 Protein

Catalog Number:PKSH030830



nucleus; where these lipids are released to a group of nuclear receptors such as peroxisome proliferator-activated receptors (PPARs).

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: techsupport@elabscience.com