

Recombinant Human CD1B/CD1A Protein (Fc Tag)

Catalog Number:PKSH030990



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	CD1;CD1A;CD1B;R1
Species	Human
Expression Host	HEK293 Cells
Sequence	Met 1-Ser 303
Accession	NP_001755
Calculated Molecular Weight	58.6 kDa
Observed molecular weight	59-69 kDa
Tag	C-hFc

Properties

Purity	> 88 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4 Normally 5 % - 8 % trehalose, mannitol and 0.01 % Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 88 % as determined by reducing SDS-PAGE.

Background

CD1B contains 1 Ig-like (immunoglobulin-like) domain and belongs to the CD1 family. CD1 family members are transmembrane glycoproteins, which are structurally related to the major histocompatibility complex (MHC) proteins and form heterodimers with beta-2-microglobulin. During protein synthesis and maturation, they bind endogenous lipids that are replaced by lipid or glycolipid antigens when the proteins are internalized and pass through endosomes, before trafficking back to the cell surface. CD1B localizes to late endosomes and lysosomes via a tyrosine-based motif in the cytoplasmic tail, and requires vesicular acidification to bind lipid antigens.. It is expressed on cortical thymocytes, epidermal Langerhans cells, dendritic cells, on certain T-cell leukemias, and in various other tissues. CD1B is an antigen-presenting protein that binds self and non-self lipid and glycolipid antigens and presents them to T-cell receptors on

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Human CD1B/CD1A Protein (Fc Tag)

Catalog Number:PKSH030990



natural killer T-cells.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: techsupport@elabscience.com