

Recombinant Human ALOX5AP/FLAP Protein (His Tag)



Catalog Number:PKSH031059

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	FLAP
Species	Human
Expression Host	Baculovirus-Insect Cells
Sequence	Met 1-Pro 161
Accession	NP_001620.2
Calculated Molecular Weight	19.0 kDa
Observed molecular weight	19 kDa
Tag	C-His

Properties

Purity	> 90 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile 20mM Tris, 150mM NaCl, 10% glycerol, pH 7.4, 0.5% triton, 5mM DTT Normally 5 % - 8 % trehalose, mannitol and 0.01 % Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 90 % as determined by reducing SDS-PAGE.

Background

Arachidonate 5-Lipoxygenase-Activating Protein (ALOX5AP), also known as FLAP, belongs to the MAPEG family. ALOX5AP/FLAP is an essential partner of 5-LO for this process. The FLAP (ALOX5AP) gene has been linked to risk for myocardial infarction, stroke and restenosis, reigniting pharmaceutical interest in this target. It had been found that ALOX5AP/FLAP is a key enzyme in leukotriene formation, in both human pulmonary microvascular endothelial cells and a transformed human brain endothelial cell line. In addition, the protein FLAP has recently been identified as an emerging target in metabolic disease. In fact, FLAP is overexpressed in the adipose tissue of patients and experimental animals with obesity.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017