

Recombinant Human PLA2G2A/PLA2B Protein (His Tag)



Catalog Number:PKSH031219

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	MOM1;PLA2;PLA2B;PLA2L;PLA2S;PLAS1;sPLA2
Species	Human
Expression Host	HEK293 Cells
Sequence	Met 1-Cys 144
Accession	NP_000291.1
Calculated Molecular Weight	15.4 kDa
Observed molecular weight	19 kDa
Tag	C-His

Properties

Purity	> 97 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4 Normally 5 % - 8 % trehalose, mannitol and 0.01 % Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 97 % as determined by reducing SDS-PAGE.

Background

Phospholipase A2, membrane associated, also known as Phosphatidylcholine 2-acylhydrolase 2A, Group IIA phospholipase A2, Non-pancreatic secretory phospholipase A2 and PLA2G2A, is a peripheral membrane protein which belongs to the phospholipase A2 family. PLA2G2A is found in many cells and also extracellularly. The membrane-bound and secreted forms of PLA2G2A are identical. PLA2G2A has been proposed to play a role in anti-bacterial defense, inflammation and eicosanoid generation, in clearance of apoptotic cells, and in the Wnt signaling pathway. PLA2G2A is thought to participate in the regulation of the phospholipid metabolism in biomembranes including eicosanoid biosynthesis. PLA2G2A catalyzes the calcium-dependent hydrolysis of the 2-acyl groups in 3-sn-phosphoglycerides. PLA2G2A might be a factor in human colorectal tumorigenesis.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017