

Recombinant Human NEGR1 Protein (His Tag)

Catalog Number:PKSH031226



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	DMML2433;IGLON4;KILON;Ntra
Species	Human
Expression Host	HEK293 Cells
Sequence	Met 1-Gly 324
Accession	NP_776169.2
Calculated Molecular Weight	32.9 kDa
Observed molecular weight	45-48 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4 Normally 5 % - 8 % trehalose, mannitol and 0.01 % Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Neuronal Growth Regulator 1, NEGR1, also known as neurotractin, or KILON, which belongs to the immunoglobulin superfamily, IgLON family. This GPI-linked cell surface glycoprotein NEGR1 is composed of three Ig-like domains and belongs to the IgLON subgroup of neural IgSF members. It is expressed in two isoforms with apparent molecular masses of 50 and 37 kD, termed L-form and S-form, respectively. NEGR1/Neurotractin participates in the regulation of neurite outgrowth in the developing brain, and is expressed on neurites of primary hippocampal neurons. Neurotractin/KILON is a trans-neural growth-promoting factor for outgrowing axons following hippocampal denervation. KILON (kindred of IgLON) and opioid-binding cell adhesion molecule belong to the IgLON subgroup of immunoglobulin superfamily together with the limbic system-associated membrane protein and neurotrimin. The alteration of modulatory function of

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Human NEGR1 Protein (His Tag)

Catalog Number:PKSH031226



KILON/NEGR1 for the number of dendritic synapses concomitant with changes in its localization and detergent solubility during neuronal culture development. In addition to its reported role in the brain, NEGR1 is also expressed in subcutaneous adipose tissue and acts as a central 'hub' in an obesity-related transcript network.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: techsupport@elabscience.com