

Recombinant Human CD8 beta/CD8B Protein (Fc Tag)

Catalog No. PKSH031296

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	T-Cell Surface Glycoprotein CD8 Beta Chain;CD8b;CD8B;CD8B1;LEU2;LY3;LYT3;P37
Species	Human
Expression Host	HEK293 Cells
Sequence	Met 1-Pro 170
Accession	P10966-1
Calculated Molecular Weight	44.0 kDa
Observed molecular weight	50 kDa
Tag	C-hFc
Bioactivity	Not validated for activity

Properties

Purity	> 93 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4 Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 93 % as determined by reducing SDS-PAGE.

Background

CD8B (CD8b molecule); also known as P37 and LEU2; contains 1 Ig-like V-type (immunoglobulin-like) domain. The

For Research Use Only

CD8 antigen is a cell surface glycoprotein found on most cytotoxic T lymphocytes that mediates efficient cell-cell interactions within the immune system. The CD8 antigen; acting as a coreceptor; and the T-cell receptor on the T lymphocyte recognize antigens displayed by an antigen presenting cell (APC) in the context of class I MHC molecules. The functional coreceptor is either a homodimer composed of two alpha chains; or a heterodimer composed of one alpha and one beta chain. Both alpha and beta chains share significant homology to immunoglobulin variable light chains. P37 gene encodes the CD8 beta chain isoforms. Multiple alternatively spliced transcript variants encoding distinct membrane associated or secreted isoforms have been described. A pseudogene; also located on chromosome 2; has been identified. CD8 is thought to play a role in the process of T-cell mediated killing.