

Recombinant Human TCN2 Protein (His Tag)

Catalog Number:PKSH031521



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Transcobalamin-II;D22S676;D22S750;TC-2;TC2;TCII
Species	Human
Expression Host	CHO Stable Cells
Sequence	Met 1-Trp 427
Accession	NP_000346.2
Calculated Molecular Weight	46.7 kDa
Observed molecular weight	44 kDa
Tag	C-His

Properties

Purity	> 90 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4 Normally 5 % - 8 % trehalose, mannitol and 0.01 % Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 90 % as determined by reducing SDS-PAGE.

Background

Transcobalamin II, also known as TCN2 and TC II, is a plasma protein that binds cobalamin (Cbl; vitamin B12) as it is absorbed in the terminal ileum and distributes to tissues. The circulating transcobalamin II-cobalamin complex binds to receptors on the plasma membrane of tissue cells and is then internalized by receptor-mediated endocytosis. Transcobalamin II is a non-glycosylated secretory protein of molecular mass 43 kDa. Its plasma membrane receptor (TC II-R) is a heavily glycosylated protein with a monomeric molecular mass of 62 kDa. Human TCN2 gene is composed of nine exons and eight introns spanning approximately 20 kb with multiple potential transcription start sites. A number of genetic abnormalities are characterized either by a failure to express TCN2 or by synthesis of an abnormal protein. The TCN2 deficiency results in cellular cobalamin deficiency, an early onset of megaloblastic anaemia, and neurological

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Human TCN2 Protein (His Tag)

Catalog Number:PKSH031521



abnormalities.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: techsupport@elabscience.com