## Recombinant Rat EphA3 Protein (Fc Tag)

Catalog Number: PKSR030159



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description	
Synonyms	EPHA3
Species	Rat
Expression Host	HEK293 Cells
Sequence	Met1-His541
Accession	EDL75897.1
Calculated Molecular Weight	85.8 kDa
Observed molecular weight	103 kDa
Tag	C-hFc
Bioactivity	Immobilized ratEFNA5-His at 10 $\mu$ g/ml (100 $\mu$ l/well) can bind rat EPHA3-Fc, The EC50 of rat EPHA3-Fc is 12.1-28.4 ng/ml.
Properties	
Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per $\mu$ g of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4 Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.
Data	
KDa M 116 66.2 45.0 35.0	

> 95 % as determined by reducing SDS-PAGE.

25.0 18.4 14.4

## **Background**

EPHA3 gene belongs to the ephrin receptor subfamily of the protein-tyrosine kinase family. EPH and EPH-related receptors have been implicated in mediating developmental events, particularly in the nervous system. The ephrin receptors are divided into 2 groups based on the similarity of their extracellular domain sequences and their affinities for binding ephrin-A and ephrin-B ligands. EPHA3 gene encodes a protein that binds ephrin-A ligands. EPHA3 is involved in the retinotectal mapping of neurons. It may also control the segregation but not the guidance of motor and sensory axons during neuromuscular circuit development.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine Toll-free: 1-888-852-8623 Tel: 1-832-243-60