## Recombinant Mouse CD39/ENTPD1 Protein (His Tag)

Catalog No. PKSM040702

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description	
Synonyms	Ectonucleoside triphosphate diphosphohydrolase 1;NTPDase 1;Ecto- ATP diphosphohydrolase 1;Ecto-ATPDase 1;Ecto-ATPase 1;Ecto- apyrase;Lymphoid cell activation antigen;CD39
Species	Mouse
Expression Host	Baculovirus-Insect Cells
Sequence	Thr 38-Ile 478
Accession	P55772
Calculated Molecular Weight	51.0 kDa
Tag	C-His
Bioactivity	Measured by its ability to hydrolyze the 5'phosphategroups from the substrate adenosine 5'triphosphate(ATP). The specific activity is > 25, 000 pmoles/min/ $\mu$ g.
Properties	
Properties Purity	> 95 % as determined by reducing SDS-PAGE.
Properties Purity Endotoxin	<ul> <li>&gt; 95 % as determined by reducing SDS-PAGE.</li> <li>&lt; 1.0 EU per μg of the protein as determined by the LAL method.</li> </ul>
Properties Purity Endotoxin Storage	<ul> <li>&gt; 95 % as determined by reducing SDS-PAGE.</li> <li>&lt; 1.0 EU per µg of the protein as determined by the LAL method.</li> <li>Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at &lt; -20°C for 3 months.</li> </ul>
Properties Purity Endotoxin Storage Shipping	<ul> <li>&gt; 95 % as determined by reducing SDS-PAGE.</li> <li>&lt; 1.0 EU per µg of the protein as determined by the LAL method.</li> <li>Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at &lt; -20°C for 3 months.</li> <li>This product is provided as lyophilized powder which is shipped with ice packs.</li> </ul>
Properties Purity Endotoxin Storage Shipping Formulation	<ul> <li>&gt; 95 % as determined by reducing SDS-PAGE.</li> <li>&lt; 1.0 EU per µg of the protein as determined by the LAL method.</li> <li>Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at &lt; -20°C for 3 months.</li> <li>This product is provided as lyophilized powder which is shipped with ice packs.</li> <li>Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 7.4, 10% glycerol Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization.</li> <li>Please refer to the specific buffer information in the printed manual.</li> </ul>
Properties Purity Endotoxin Storage Shipping Formulation Reconstitution	<ul> <li>&gt; 95 % as determined by reducing SDS-PAGE.</li> <li>&lt; 1.0 EU per µg of the protein as determined by the LAL method.</li> <li>Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at &lt; -20°C for 3 months.</li> <li>This product is provided as lyophilized powder which is shipped with ice packs.</li> <li>Lyophilized from sterile 20mM Tris, 500mM NaCl, pH 7.4, 10% glycerol Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization.</li> <li>Please refer to the specific buffer information in the printed manual.</li> </ul>

Data



> 95 % as determined by reducing SDS-PAGE.

## Background

## **For Research Use Only**

## **Elabscience**®

CD39, also known as ENTPD1, belongs to the GDA1/CD39 NTPase family. It is expressed primarily on activated lymphoid cells and can also be detected in endothelial tissues. The vascular isoform and the placental isoform II are present in both placenta and umbilical vein, whereas placental isoform I is present in placenta only. CD39 can hydrolyze both nucleoside triphosphates and diphosphates. It is the dominant ecto nucleotidase of vascular and placental trophoblastic tissues and appears to modulate the functional expression of type 2 purinergic (P2) G protein coupled receptors (GPCRs). CD39 transgenic mice exhibit impaired platelet aggregation, prolonged bleeding times, and resistance to systemic thromboembolism. There is a correlation between ATP hydrolysis and triglycerides in patients with chronic heart disease, suggesting a relationship between ATP diphosphohydrolase and thrombogenesis. In the nervous system, CD39 could hydrolyze ATP and other nucleotides to regulate purinergic neurotransmission.

**For Research Use Only**