

Recombinant DENV (type 1, strain US/Hawaii/1944) DENV-NS1 Protein (His Tag)

DIA • AN[®]
by Elabscience

Catalog Number:PKSV030122

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	DENV-NS1 Protein;DENV;Dengue virusNS1 protein
Species	DENV
Expression Host	HEK293 Cells
Sequence	Asp776-Ala1127
Accession	ACF49259.1
Calculated Molecular Weight	42.2 kDa
Tag	N-His

Properties

Purity	> 90 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 90 % as determined by reducing SDS-PAGE.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017