

Recombinant DENV (type 2, strain New Guinea C/PUO-218 hybrid) E / Envelope Protein (His Tag)



Catalog Number:PKSV030128

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	E Protein;DENV;Dengue virus
Species	DENV
Expression Host	Baculovirus-Insect Cells
Sequence	Met281-Gly679
Accession	AAC59274.1
Calculated Molecular Weight	45.7 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per μg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C . Reconstituted protein solution can be stored at $4-8^{\circ}\text{C}$ for 2-7 days. Aliquots of reconstituted samples are stable at $< -20^{\circ}\text{C}$ for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile 20 mM Tris, 300 mM NaCl, pH 8.5, 10 % glycerol. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017